



About your hearing aids



Hearing aid controls

Your right hearing aid has a red marker.

Your left hearing aid has a blue marker.



Close the battery drawer to turn your hearing aid on.




Open the battery drawer to turn your hearing aid off.




Press the top of the button to turn the volume up.

Press the bottom of the button to turn the volume down.



Batteries

Batteries last 7-14 days.
Change the battery when it runs out



The hearing aid will play a tune
when the battery is low




Take the old battery out



Take the sticker off the new battery



Put the new battery in your
hearing aid



How to get new batteries

Batteries come in packets like this



You can get new batteries from:

You must take your battery card



Putting your hearing aid in

Put the earmould into your ear



Put the hearing aid behind your ear



You are now ready



What to expect


Your hearing aid might sound strange at the start




You will soon get used to the sound



Wear it all day



Take it out at night



The Loop System


You can use the loop when you see this picture



Here are some places that might have the loop



Press and hold the button to connect to the loop



Press and hold the button to go back to the normal setting



Looking after your hearing aids


Keep your earmould clean



Visit us every 6 months to change the tubes



Tell us if there is a problem



How to contact us:

Phone: 01392 402223



Email: cse.audiology@nhs.net



Website: www.chimehealth.co.uk



Document control

About Stretta Ms	Latest review date: 07-2024	Version number: 1.1
Status: Approved	Next review: 07-2025	Last page