




# About your hearing aids



# Hearing aid controls

Your left hearing aid has a blue marker.

Your right hearing aid has a red marker.



Close the battery drawer to turn your hearing aid on.




Open the battery drawer to turn your hearing aid off.




Press the top of the button to turn the volume up.

Press the bottom of the button to turn the volume down.



## Batteries

Batteries last 7-14 days.  
Change the battery when it runs out



The hearing aid will play a tune when the battery is low



Take the old battery out



Take the sticker off the new battery



Put the new battery in your hearing aid



## How to get new batteries

Batteries come in packets like this



You can get new batteries from:

You must take your battery card



**My hearing aid battery card**

Name: \_\_\_\_\_

Postcode: \_\_\_\_\_

## Putting your hearing aid in

Put the hearing aid behind your ear



Push the dome into your ear



You are now ready



## What to expect


Your hearing aid might sound strange at the start




You will soon get used to the sound



Wear it all day



Take it out at night




## Looking after your hearing aids


Keep your hearing aid clean



Visit us every 6 months to change the domes



Tell us if there is a problem



## How to contact us:


Phone: 01392 402223



Email: [cse.audiology@nhs.net](mailto:cse.audiology@nhs.net)



Website: [www.chimehealth.co.uk](http://www.chimehealth.co.uk)



### Document control

About Stretta RICs	Latest review date: 07-2024	Version number: 1.1
Status: Approved	Next review: 07-2025	Last page